

## January Release Notes – 01/08/2018

### Contents

Feature Enhancements .....	1
Readmission of deceased patients .....	1
Advanced Search.....	1
Feature Updates .....	2
Manage Directory .....	2

### **Feature Enhancements**

#### **Readmission of deceased patients**

GEHRIMED now enables the readmission of patients with an erroneous 'Death' status. When creating an encounter for a deceased patient with an encounter date after the date they were marked deceased, the following prompt will appear:



The screenshot shows a web interface titled "Create Encounter". At the top right, there are help and close icons. Below the title bar is a red error message box that says "Please correct the errors below and try again." Below this, there is a form with a "Patient" label and a text input field containing "EXAMPLE PATIENT 1". To the right of the input field is an "Advanced Search" button. Below the input field, there is a red error message: "The patient record indicates the patient is deceased as of 12/05/2017." Below this message is a checkbox with the text "This patient is not deceased, and I confirm resetting the patient's status to active." Below the checkbox is a blue link that says "Create New Patient". At the bottom of the form, there is a search bar with a magnifying glass icon and a "16" next to it.


By selecting the checkbox, you confirm that this patient is not actually deceased, and by creating the encounter you are changing their status back to 'Active'.

#### **Advanced Search**

The GEHRIMED advanced search has been updated with a number of optimizations to improve search performance:

#### **Pagination**

Search results will now be divided into pages, defaulting to 10 results per page:



You may select the number of results you want displayed per page, up to the maximum of 500 (GEHRIMED limits search results to 500 displayed entries).


## Printing

The Advanced Search Print function will now print only the results displayed on the page. To print all displayed entries, you must select to display 500 per page. To access results over 500, you must choose to export the search results. Export functionality has not been changed, and will still export up to 10,000 results.

## Feature Updates

### Manage Directory

The Manage Directory listings will now be tiered to display Last Name, First Name. If the Last and First Names are unavailable, the Practice Name will be displayed instead:



~end of document~