

Release Notes – June 2, 2021

The GEHRIMED iOS application will be updated for functionality improvements.

NOTE: This update applies only to iPad users. Microsoft Windows users are not affected.

Issues Addressed

- iPad issues relating to voice recognition software
 - Crashes after startup, when entering diagnoses or when switching between fields
 - Microphone button functionality improvements
- The newest version of GEHRIMED is 1.3.24.10
- GEHRIMED recommends uninstalling current version of iOS and reinstalling iOS 14.5.1 or higher for best performance
- Additional user resources for speech-to-text functionality will be available in Zendesk
- **NOTE:** To view your currently installed version of iOS or GEHRIMED, open the application and check the bottom right corner of window:



~end of document~